

Site Assessment Form

Gathering Information and Specifying System

Measuring a room for a Ceiling Lift

Measure and specifying a ceiling lift system is not a difficult task, and can be accomplished by one person. The following items below are needed to measure a room for a lift system:

- 25' Tape Measure
- Site Assessment Form
- Flashlight
- Small 4' Step Ladder
- Digital Camera and/or Video Camera

Once familiar with how a system is installed, the time required to measure a system will be reduced and the information gathering will become easier. As a general rule, more information is better, and digital pictures are a great help.

The type of mounting method (wall support vs ceiling mount) will alter the information required, and understanding the differences can allow for some information to be omitted. If the installation method is not known at the time of site, all information should be gathered.

Territory Contact Information

Dealer Name:	Phone:
Dealer Rep Name:	

Project Information

Project Address:		
City:	State:	Zip Code:

Ceiling Lift Information


Fixed Lift		Portable Lift	
Load Capacity	Traverse/Charging Type	Load Capacity	Portable Accessories:
<input type="checkbox"/> 300 lbs.	<input type="checkbox"/> Manual	<input type="checkbox"/> 300 lbs.	<input type="checkbox"/> 24" Reacher Arm
<input type="checkbox"/> 450 lbs.	<input type="checkbox"/> Power	<input type="checkbox"/> 440 lbs.	<input type="checkbox"/> 36" Reacher Arm
<input type="checkbox"/> 625 lbs.	<input type="checkbox"/> Power Return to Charge	<input type="checkbox"/> 600 lbs.	Transport Cart: Y or N
<input type="checkbox"/> 800 lbs.	<input type="checkbox"/> Standard (Docking Station)		
<input type="checkbox"/> 1000 lbs.	<input type="checkbox"/> Omni (Constant Charge)		

Sling Information

Client Height:	_____ feet, _____ inches
Client Weight:	
Head Support:	<input type="checkbox"/> Yes <input type="checkbox"/> No

PRISM MEDICAL

Please provide a sketch of the ceiling track layout with dimensions. Include walls, doorways, and ceiling obstructions if applicable.



Please indicate joist direction:

